
Low-cost biosensors with printed interdigitated graphene electrodes

Thomas Velten, Thorsten Knoll, Axel Brenner, Anke Schultz (*Fraunhofer IBMT*)

Joachim Wiest (*cellasys GmbH - R&D*)

Renate Warmers, Gerald Jenke, Anna Zumbülte (*Saueressig GmbH + Co. KG*)

Andreas Urban (*AiCuris Anti-infective Cures GmbH*)


Kris Seunarine (*Haydale Ltd.*)

Outline

- Motivation
 - What do we want to measure
 - Drawback of the state of the art
- BIOGRAPHY approach
- Results
 - Printing
 - Biology
- Summary and outlook


What do we want to measure

- Sensor for toxicity testing
 - Reaction of living cells to exposure to substances
- Sensor for validation of antiviral substances
 - Cells insulate electrodes
 - Addition of viruses → viruses destroy cells (cytopathic effect ,CPE')
 - Addition of viruses and antiviral substances → CPE is inhibited




Drawback of commercial systems

- High costs due to
 - Production process for printed circuit boards
 - Biocompatible metals (gold, platinum)





→ Cost effective sensor production method is required


BIOGRAPHY approach



+



- Electrically conductive
- Biocompatible
- Improves cell adhesion




Roll-to-roll gravure printing machine R²R-300



Roll-to-roll gravure printing machine R²R-300


- Two printing units: each with an integrated NIR drying unit
- Power of corona unit: Up to 600 W
- Printing speed: Up to 40 m/min
- Web width: 300 mm
- Width of printable area: 250 mm
- Doctor blade pressure and angle are independently adjustable



Gravure printing cylinders


Cylinder for graphene printing

- Masking and chemical etching
- → 50 μm deep cells, 80 l/cm; high ink transfer (low electrical resistance)



Cylinder for protein printing

- Direct engraving in copper with ultrashort pulse laser
- → microstructures < 10 μm
- Alternative: 20 μm deep cells; 120 l/cm



Directly structured dots
with diameter 8 μm

Graphene ink

- Basis: screen printing graphene ink (HDPlas IGSC02002)
 - Surface-functionalised Graphene Nanoplatelet (GNP) material (using Haydale's patented HDPlas® plasma-based technology)
 - Viscosity: $\sim 5.5 \text{ Pa}\cdot\text{s}$
 - Thickness before / after drying: $\sim 13 \mu\text{m} / 7 \mu\text{m}$
 - Sheet resistivity: $\sim 12 \Omega/\text{sq.}$ (normalised to $25 \mu\text{m}$)




- Modifications
 - Solvent content
 - Resin content
 - Milling procedure

Protein ink

- Composition
 - Matrigel, dissolved in DMEM/F12
 - Ammonium acetate buffer
 - Carboxymethylcellulose sodium salt, dissolved in DI water
 - X

Roll-to-roll gravure printing with graphene ink



- Printed graphene patterns with test cylinder



Electrical characterization


- Resistivity measurements
- Line width: ~1.05 mm, length: 10 mm
- Angles: 45° and 135°, cell depth: 60 μm
- Graphene thickness: ~ 4 μm

$$R_{sq.} = 10 - 20 \Omega/sq.$$



Cell adhesion and proliferation on graphene electrodes

- FDA (fluorescein diacetate) stained TZM-bl cells on graphene after 48 h incubation



Cytotoxicity testing of graphene ink (ISO 10993)

- Cytotoxicity of PET foil as well as of graphene ink printed PET foil was analysed in more than ten different cell lines
- Cells are from
 - different species (human, mouse, rat)
 - different organs and tissues (liver, kidney, heart, nerve system, skin and blood)
- No evidence that PET foil or graphene ink has any cytotoxic potential

Final printing parameters (graphene ink)

- PET foil: thickness 50 μm , surface energy < 40 mN/m
- Surface energy of graphene ink ~ 45 mN/m
- Corona activation (600 W) \rightarrow surface energy ~ 60 mN/m


- Speed: 40 m/min
- Doctor blade contact pressure: 4.0 bar
- Cylinder contact pressure: 3.5 bar
- NIR drying power: 90% (of 3.5 kW)

Final printing parameters (protein ink)


- Speed: 40 m/min
- Doctor blade contact pressure: 3.0 bar
- Cylinder contact pressure: 2.5 bar
- NIR drying power: 60% (of 3 kW)

Printed structures


120 l/cm
30 μ m



A (52/52)






B (100/100)







B1 (100/200)


80 l/cm
50 μ m




96-well plate set-up with printed circuit adapter




96-well plate set-up with printed circuit adapter



Scale bar: 3.5 mm



Drift after adding cell culture medium at 37°C



Stable after < 10 hours

Summary

- We have successfully developed
 - Graphene ink for gravure printing
 - Protein ink for gravure printing
 - Two-colour rotogravure printing system
 - Processes for printing graphene electrodes and for printing proteins
- Printed structures are biocompatible
- Printed proteins improve cell adhesion/proliferation
- Printed sensors are useable with the IMOLA - IVD system of cellasys - R&D
- **Next steps**
 - Validation of biosensors for the two application fields

Acknowledgements

SPONSORED BY THE



Federal Ministry
of Education
and Research

MANAGED BY



PTKA

Project Management Agency Karlsruhe

Karlsruhe Institute of Technology

Innovate UK
Technology Strategy Board

This research and development project is funded by the German Federal Ministry of Education and Research (BMBF) within the Framework Concept "Research for Tomorrow's Production" (funding number 02PN2240) and managed by the Project Management Agency Karlsruhe (PTKA). The author is responsible for the contents of this publication.

Furthermore, this project is co-funded by Innovate UK. Innovate UK is the UK's innovation agency. It works with people, companies and partner organisations to find and drive the science and technology innovations that will grow the UK economy.

Thank you

Visit us at the booth of Saueressig GmbH + Co.KG

(Hall B0 Booth 313)

BIOGRAPHY consortium

M-Era.net project with partners from Germany and UK


Tasks	Name
R2R gravure printing process (graphene and proteins), biocompatibility tests	Fraunhofer IBMT
R2R two-colour printing machine, laser-based method for micro patterning of gravure printing cylinders	SAUERESSIG GmbH + Co. KG
Biocompatible graphene ink for gravure printing	Haydale Ltd.
Application: antiviral drug validation	AiCuris GmbH & Co. KG
Application: repeated dose toxicity	cellasys GmbH R&D



Investigation of impact of the graphene ink on the inhibitory effect of antiviral drugs

- Serial dilutions of two different antivirals were pre-incubated with graphene printed PET foils before a regular antiviral assay was performed
- Antiviral drugs used:
 - BAY 41-4109 (inhibitor for HBV)
 - Efavirenz (inhibitor for HIV)
- There was no evidence that graphene ink has any influence on the antiviral activity of BAY 41-4109 and Evavirenz

Repeated annealing in air at 50 °C



Mechanical wear and impedance

- 5x insertion of sensors into connector to test dependence of impedance values on mechanical wear

